



COLTENE

restorative product guide



ONE COAT 7 UNIVERSAL

SINGLE COMPONENT LIGHT CURED UNIVERSAL ADHESIVE

ONE COAT 7 UNIVERSAL is a light-cured, one-component bonding agent used with the self-etching, selectively etching or total etch technique for adhesive restorations.



- Offers the ability to Self-Etch on the preparation
- Contains the monomer MDP for optimal adhesion to dentin, enamel, and metals
- Can be used in combination with light cure resin composite or dual cure self-adhesive resin cement, SoloCem
- 5 ml bottle design offers ~300 applications

ONE COAT 7 UNIVERSAL Order Information

60019540 ONE COAT 7 UNIVERSAL Single Dose
50 × 0.1 ml ONE COAT 7 UNIVERSAL Single Dose
50 × Microbrush, black

60019539 ONE COAT 7 UNIVERSAL Refill
1 × 5 ml ONE COAT 7 UNIVERSAL refill

60019538 ONE COAT 7 UNIVERSAL Intro Kit
1 × 5 ml One Coat 7 Universal
1 × 2.5 ml Etchant Gel S
1 × 20 pcs. Application tips
1 × 50 pcs. Brushes, black
1 × Mixing well
1 × Cover for mixing well



Solo Cem[®]

SELF-ADHESIVE DUAL-CURING RESIN CEMENT

SoloCem is a universal cement that is safe for luting with or without additional bonding. An all-in-one cement for use in almost all indirect luting indications, it is easy to use, provides long-lasting restorations and requires less inventory.



- Contains the monomers 4-META and MDP for optimal adhesion to dentin, enamel, and metals
- Can be used as is or in combination with ONE COAT 7 UNIVERSAL for improved bond strength
- In direct restorations:
 - ◆ Retentive Preps / PFM & Zirconia Restoration – SoloCem
 - ◆ Smooth Preps / All Ceramic & Lithium Disilicate Restoration – ONE COAT 7 UNIVERSAL on Prep + SoloCem
- Easy clean up, low film thickness, and high radiopacity of 4mm Al/mm
- Contains zinc oxide



SoloCem Order Information

- | | |
|------------------|--|
| 6001 3986 | SoloCem Intro Kit |
| 1 × 5ml | SoloCem Dentin |
| 1 × 5ml | SoloCem Translucent |
| 10 × | Mixing Tip Brown Short Super Fine |
| 10 × | Mixing Tip Brown Short Fine |
| 60014062 | SoloCem Dentin Refill |
| 1 × | SoloCem Dentin, 5 ml |
| 5 × | Mixing Tip Brown Short Super Fine |
| 5 × | Mixing Tip Brown Short Fine |
| 60014061 | SoloCem White Opaque Refill |
| 1 × | SoloCem White Opaque, 5 ml |
| 5 × | Mixing Tip Brown Short Super Fine |
| 5 × | Mixing Tip Brown Short Fine |
| 60019840 | SoloCem Trans Refill |
| 1 × | SoloCem Trans, 5 ml |
| 5 × | Mixing Tip Brown Short Super Fine |
| 5 × | Mixing Tip Brown Short Fine |
| 6747 | Mixing Tip Brown Short Super Fine |
| 40 × | Mixing Tip, Ø 1 mm |
| 6759 | Mixing Tip Brown Short Fine |
| 40 × | Mixing Tip, Ø 1.8mm |



BRILLIANT EverGlow[®]

UNIVERSAL SUBMICRON HYBRID COMPOSITE

Invented by COLTENE, Duo Shades simplify shade management by combining two VITA* shades into a single shade (for example, A1/B1 or A2/B2). This is made possible by the optical properties of Duo Shades, which allow them to blend in harmoniously with adjacent teeth.

- Exceptional polishability and long-lasting brilliance
- Aesthetic, single-shade restorations
- Ideal handling
- Good wettability
- Contains zinc oxide

“I was surprised how much I ended up liking this composite. It had great color matching, excellent shades available, and a great polish.”

Nathan Vassiliades, DMD

* Not a registered trademark of VITA



BRILLIANT EverGlow® Flow

SUBMICRON HYBRID FLOWABLE COMPOSITE

BRILLIANT EverGlow® Flow is a submicron hybrid flowable composite that offers both stability and flowability in one versatile material. User friendly, aesthetic, and appropriate for use in a wide range of indications, it is the perfect complement to BRILLIANT EverGlow®.

- Suitable for initial layers in all cavity classes
- Convenient Duo Shades
- Good wettability for extensive coverage
- Two application needle sizes for multiple indications
- Contains zinc oxide



/// The choice of tips made BRILLIANT EverGlow Flow versatile for filling smaller hard-to-reach areas on roots and for use as a liner/base. ///

Jeff Green, DDS

Brilliant EverGlow Order Information

Refills Tips (20 × 0.2 g)

60019746	Opaque Bleach (OBL)
60019747	Opaque A1 (OA1)
60019748	Opaque A3 (OA3)
60019706	Bleach (BL)
60019707	A1 / B1
60019708	A2 / B2
60019709	A3 / D3
60019710	A3.5 / B3
60019716	C2 / C3
60019717	A4 / C4
60019718	Translucent (Trans)
60019719	Bleach Translucent (BL Trans)

Refills Syringes (1 × 3 g)

60019743	Opaque Bleach (OBL)
60019744	Opaque A1 (OA1)
60019745	Opaque A3 (OA3)
60019697	Bleach (BL)
60019698	A1 / B1
60019699	A2 / B2
60019700	A3 / D3
60019701	A3.5 / B3
60019702	C2 / C3
60019703	A4 / C4
60019704	Translucent (Trans)
60019705	Bleach Translucent (BL Trans)

Accessories

60021919	BRILLIANT EverGlow Shade Guide
60022485	Composite Tip Dispenser

Brilliant EverGlow Flow Order Information

Refills (1 × 2 g + 10 application needles)

60019751	Opaque A2 (OA2)
60019752	Bleach (BL)
60019753	A1 / B1
60019754	A2 / B2
60019755	A3 / D3
60019756	A3.5 / B3
60019757	A4 / C4
60019758	Translucent (Trans)

Accessories

8469	20 Application needles Standard (Ø 0.8 mm)
60022137	20 Application needles Fine (Ø 0.4 mm)






BRILLIANT Crios

REINFORCED COMPOSITE BLOCK FOR PERMANENT RESTORATION

BRILLIANT Crios is the ideal choice for single tooth restorations, both in the anterior and the posterior region.

- Reliable, aesthetic restorations without a separate firing process
- Indicated for inlays, onlays, crowns and veneers
- Shock absorbing effect



Suitable for inlays, onlays, crowns and veneers. Available for CEREC/InLab devices.

Brilliant Crios Order Information

Low Translucent 12

Bleach	60019998	BL LT 12
A1	60019999	A1 LT 12
A2	60020000	A2 LT 12
A3	60020001	A3 LT 12
A3.5	60020002	A3.5 LT 12
B1	60020003	B1 LT 12
B2	60020004	B2 LT 12
B3	60020005	B3 LT 12
C2	60020006	C2 LT 12

Low Translucent 14

Bleach	60019985	BL LT 14
A1	60019986	A1 LT 14
A2	60019987	A2 LT 14
A3	60019988	A3 LT 14
A3.5	60019989	A3.5 LT 14
B1	60019990	B1 LT 14
B2	60019991	B2 LT 14
B3	60019992	B3 LT 14
C2	60019993	C2 LT 14

High Translucent 12

A1	60020007	A1 HT 12
A2	60020008	A2 HT 12
A3	60020009	A3 HT 12
B1	60020010	B1 HT 12

High Translucent 14

A1	60019994	A1 HT 14
A2	60019995	A2 HT 14
A3	60019996	A3 HT 14
B1	60019997	B1 HT 14



Fill-Up!™

DUAL CURING BULK COMPOSITE

This medium viscosity, dual-cure bulk composite has been developed for those indications that demand fast, aesthetic tooth-colored restorations.

- Optimal depth polymerization with minimal shrinkage due to dual-cure system
- Single layer technique even in very deep preparations
- Optimal marginal sealing with the matched ParaBond bonding agent
- Universal shade in a convenient automix syringe for efficient placement
- Contains zinc oxide

FAST RESULTS FOR CLASS I AND CLASS II CAVITIES



Apply in a single step



Light cure for 5 seconds or wait for the chemical setting



Adapt and polish the filling



Final result

It was very easy to place in large restorations due to good flow and the delivery system. It made it a lot easier in hard-to-reach posterior areas.

Dr. Milton Ruiz

Fill-Up! Order Information

60019341 Fill-Up! Refill
1 × 4.5 g Fill-Up!™ Universal
3 Mixing Tip Short Super Fine
4 Mixing Tip Short Fine

60019342 Fill-Up!™ Intro Kit
2 × 4.5 g Fill-Up!™ Universal
1 × 2.5 ml Etchant Gel S
2 × 3 ml ParaBond Adhesive A & B
6 × Mixing Tip Short Super Fine
8 × Mixing Tip Short Fine



Optimal marginal sealing with ParaBond bonding agent



ParaCore

DUAL CURING CORE AND RESIN CEMENT

ParaCore is a glass-reinforced composite, which is radiopaque and available in three shades—dentin, white and translucent. This integrated bond and cement system complements a wide array of endodontic post systems for clinically sound, long-lasting post and core restorations.

- Three Shades: dentin, white and translucent
- Optimal monoblock bond interface between the dentin-post-crown
- Superior bond strength and durability
- Stable consistency
- Contains fluoride
- Radiopaque

ParaCore material matches closely in shade to dentin and is very strong. It also saves time when used as a cement and core build-up material. The entire system works in harmony to reinforce endodontically treated teeth that require a foundation build-up.

Dr. James Braun

ParaCore Order Information

5885 ParaCore 5 ml Intro Kit

- 1 × 5 ml ParaCore Dentin
- 1 × 5 ml ParaCore White
- 1 × 5 ml ParaCore Dentin SLOW
- 1 × 5 ml ParaCore Trans SLOW
- 3 × 3 ml ParaBond Non-Rinse Conditioner + Adhesive A & B
- 20 × Mixing Tip Brown Short Super Fine

5886 ParaCore White Refill

- 2 × 5 ml ParaCore White
- 20 × Mixing Tip Brown Short Super Fine

5887 ParaCore Dentin Refill

- 2 × 5 ml ParaCore Dentin
- 20 × Mixing Tip Brown Short Super Fine

60011391 ParaCore Dentin SLOW Refill

- 2 × 5 ml ParaCore Dentin SLOW
- 20 × Mixing Tip Brown Short Super Fine

60011392 ParaCore Trans SLOW Refill

- 2 × 5 ml ParaCore Trans SLOW
- 20 × Mixing Tip Brown Short Super Fine

5853 ParaCore 25 ml Dentin Refill

- 1 × 25 ml ParaCore Dentin
- 20 × Mixing Tip Yellow / Oral Tip Yellow

5854 ParaCore 25 ml White Refill

- 1 × 25 ml ParaCore White
- 20 × Mixing Tip Yellow / Oral Tip Yellow

7486 ParaBond Kit

- 1 × 3 ml ParaBond Non-Rinse Conditioner
- 1 × 3 ml ParaBond Adhesive A
- 1 × 3 ml ParaBond Adhesive B

Accessories

- 4470 Coltène Dispenser 25 ml



Alpen® ShapeGuard

POLISHERS

ShapeGuard facilitates intraoral finishing and polishing of the following areas: interproximal, occlusal, incisal, labial, buccal, and lingual

- Morphology of tooth is preserved
- Fast, easy and intuitive finishing & polishing
- Optimized polishing systems for composite, ceramic, and zirconia
- Outstanding high luster polishing results

Alpen ShapeGuard Order Information

Composite Plus



Dia. in 1/10 mm	110	140	110	140
Head L in mm	1.60	1.60	1.60	1.60
Polisher type	Spiral Wheel			

RA/ 204

Order No.	60021904	60021905	60021906	60021907
Content	5	5	5	5

60021910 ShapeGuard Composite Plus Assorted Pack

Content	1	1	1	1
---------	---	---	---	---

Ceramic Plus



Dia. in 1/10 mm	140	140	140
Head L in mm	1.60	1.60	1.60
Polisher type	Spiral Wheel		

RA/ 204

Order No.	60021911	60021912	60021913
Content	5	5	5

60021914 ShapeGuard Ceramic Plus Assorted Pack

Content	1	1	1
---------	---	---	---

Zirconia Plus



Dia. in 1/10 mm	140	140
Head L in mm	1.60	1.60
Polisher type	Spiral Wheel	

RA/ 204

Order No.	60022205	60022206
Content	5	5

60022207 ShapeGuard Zirconia Plus Assorted Pack

Content	2	2
---------	---	---

S.P.E.C. 3

LED CURING LIGHT

Performance without compromise.

- Cures 2 mm of composite in 5 seconds (Standard Mode)
- Cures 2 mm of composite in 1 second (3K Mode)
- Sealed construction of anodized aluminum casing ensures quick and easy disinfection between patients
- Can perform 300 ten-second cures between charges (standard mode)
- Multiple curing modes, including ortho mode, for maximum functionality

S.P.E.C. 3 Order Information

60013941	S.P.E.C. 3 LED Unit USA Plug
60013953	S.P.E.C. 3 Charger Base
60013955	Power Supply
60013950	Clear light guide Turbo-Tip 8 mm
60013951	Clear light guide 11 mm
60019326	Black light guide Turbo-Tip 8 mm
60019327	Black light guide 11 mm
60013949	Barrier Sleeves, Package of 100
60013952	Battery Module



Coltolux LED

LED CURING LIGHT

Compact and efficient; the workhorse of your operatory

- Slim design allows posterior reach with better visibility
- Cures 2 mm of composite in 10 seconds
- 20-second curing cycle with an audible signal every 10 seconds
- Dome-shaped lens concentrates light output on the restoration
- Advanced reflector design minimizes light diffusion
- Minimum of 300, 20-second curing cycles with a fully charged battery

Coltolux LED Order Information

C7970100115 Coltolux LED Curing Light

1x Handpiece
1x Charger base
1x Power supply
3x Light shields
5x Light lenses
100x Barrier sleeves
1x Curing test fixture

C7957 1000 Barrier Sleeves



Coltène/Whaledent Inc.
235 Ascot Parkway
Cuyahoga Falls, Ohio 44223 / USA
T +1 330 916 8800
F +1 330 645 8704
info.us@coltene.com

Coltène/Whaledent AG
Feldwiesenstrasse 20
9450 Altstätten / Switzerland
Tel. +41 (0)71 757 53 00
Fax +41 (0)71 757 53 01
info.ch@coltene.com

Coltène/Whaledent GmbH + Co. KG
Raiffeisenstrasse 30
89129 Langenau / Germany
Tel. +49 (0)7345 805 0
Fax +49 (0)7345 805 201
info.de@coltene.com